
Freedom60®

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INTRODUCTION

The Freedom60 Infusion System consists of the following components* manufactured by KORU Medical Systems, Inc. (KORU). See **page 5** for a complete system diagram and more details.

- **Freedom60® Infusion Pump** (referred to as “Freedom60” or “pump”)
- **Freedom60® Prefilled Syringe Adapter** (referred to as “adapter”)
- **Precision™ Flow Rate Tubing** (referred to as “Precision tubing” or “tubing”)
- **Precision™ Flow Rate Controller**
- **High-Flo™ SubQ Needle Set** (referred to as “High-Flo needle set” or “needle set”)
- **Precision™ Flow Rate Controller** (referred to as “flow controller”)

**Product availability and regulatory status may vary by geographic region. Contact your healthcare provider for more information on products available in your region.*

This guide is intended for patients and/or caregivers. The information in this guide is consistent with HCP-348406 rev A. Healthcare providers should refer to the Healthcare Provider Guide for additional detail on the device.

 **CAUTION:** Patients and their caregivers must be trained by their qualified healthcare provider prior to self-administration. Patients are advised to contact their healthcare provider for all questions related to their treatment.

INTENDED PURPOSE

The Freedom60® Infusion System is intended to subcutaneously administer prescribed medicinal products to the body for therapeutic purposes.

INDICATIONS FOR USE

The Freedom60® Infusion System is indicated for adults, adolescents, and pediatrics who are prescribed infusion therapy with subcutaneously administered immunoglobulins per the infusion rates and volume per infusion site specified in the approved medicinal product labeling.

The Freedom60 Infusion System is indicated for use with the following syringes:

- Becton Dickinson & Co. BD® Plastipak™ Luer-Lok® 50 ml
- Medline 60 ml Syringe (Item Number: SKU SYR160010)
- Hizentra® 50 ml Single-Use Prefilled Syringe (Item Number: 44206-0455-25)

CONTRAINDICATIONS

The Freedom60® Infusion System is contraindicated for the delivery of blood, insulin, and medicinal products requiring greater accuracy of delivery (i.e., critical or life-sustaining medicinal products).

CLINICAL BENEFIT

The Freedom Infusion System enables a patient to self-administer subcutaneous infusion of a prescribed medicinal product regardless of setting. It provides convenience, promotes independence, and enables mobility for the patient. The Freedom Infusion System reduces the occurrence of common local site reactions such as swelling, redness, or soreness by applying constant pressure to the syringe plunger delivering the medicinal product. Under constant pressure, the flow rate slows down when increased resistance (back pressure) is encountered at the infusion site.

MRI SAFETY INFORMATION

The Freedom60 Infusion System is MR-unsafe. Do not use the Freedom60 Infusion System in or around an MRI machine.

TECHNICAL SPECIFICATIONS

Testing was performed in a controlled test lab environment at a temperature of 20-25°C and atmospheric pressure of 1,01 bar ($\pm 0,09$). System flow rate accuracy ($\pm 15\%$) was verified with a controlled medicinal product (7,4 cP PEG Surrogate) between 20-25°C.

Environmental temperature and medicinal product temperature will significantly affect the flow rate. For maximum accuracy, perform infusions at a temperature of 20-25°C. Ambient humidity and environmental pressure are not expected to affect the flow rate.

To obtain maximum accuracy of the pump, position the height of the pump within $\pm 7,6$ cm (3") of the infusion site, whether infusing in a stationary position or in motion. The following table indicates the potential variation in flow rate at different heights.

Vertical Height (inches)	% Variation From Target Flow Rate
$\pm 7,62$ cm (3 inches) from infusion site	Equivalent to Level
$\pm 15,24$ cm (6 inches) from infusion site	up to $\pm 1,2\%$ from target flow rate
$\pm 30,48$ cm (12 inches) from infusion site	up to $\pm 2,4\%$ from target flow rate
$\pm 60,96$ cm (24 inches) from infusion site	up to $\pm 4,8\%$ from target flow rate

WARNINGS

- Patient tolerance may vary. Potential side effects may include local site reactions such as swelling, redness, or soreness. If you are experiencing pain or discomfort, contact your healthcare provider. Refer to the medicinal product prescribing information for potential medicinal product reactions and side effects (e.g., paresthesia, headache, fever, dizziness, nausea, weakness, lethargy, etc.).
- Flow rates can be affected by multiple factors such as temperature, patient conditions, height differences between the system and infusion site, and variations in medicinal product viscosity. Refer to the device technical specifications for information on device performance.
- Do not connect the High-Flo needle set directly to the syringe. Connecting the needle set directly without the tubing/luer disc may cause the syringe to eject from the Freedom60 Infusion Pump and may cause internal damage to the pump.
- Do not re-sterilize or reuse the Precision tubing or High-Flo needle set. Re-using a single-use device may result in infection or malfunction of the pump.
- The black tab that pushes on the syringe plunger operates under high force. Do not place fingers on the black tab or inside the syringe shield at any time. Do not attempt to interfere with the movement of the black tab at any time.
- Do not attempt to remove the syringe or disconnect the tubing set without first turning the pump to the OFF position and fully winding the large knob clockwise until the black tab has reached the end of its track.
- The Freedom60 Infusion Pump does not have an alarm. No alarm will sound if an interruption to flow occurs. There is no display of infusion status.
- If the Freedom60 Infusion Pump is submerged in any fluid, discontinue use and call or contact your healthcare provider for a replacement.
- Do not use the Freedom60 Infusion System for blood transfusions or with medicinal product where delay or under-infusion could result in serious injury.
- The Freedom60 Infusion System does not contain a particulate filter. Do not use the Freedom60 Infusion System with medicinal products that require use of a particulate filter.

- Do not sterilize the infusion pump. Always follow the cleaning and disinfection instructions on **page 14**.
- The Freedom60 Infusion System is MR-unsafe. Do not use the Freedom60 Infusion System in or around an MRI Machine.



PRECAUTIONS

- The Freedom60 Infusion System is for prescription use only. Use the Freedom60 Infusion System only for the patient for whom the device is prescribed and only for its intended use. Ensure that the patient and/or caregiver is trained by a qualified healthcare professional prior to use.
- Use only Freedom Infusion System components manufactured by KORU Medical Systems, Inc. Use of other products may result in unknown flow rates.
- Use only compatible syringes with the Freedom60 Infusion Pump. Using a different syringe may result in unknown flow rates.
- Before use, carefully inspect the Precision tubing and High-Flo needle set packaging. Do not use the set if the package is opened.
- Inspect tubing and needle sets for damage. If damaged, do not use. Replace and contact your healthcare provider.
- The slide clamp included on the Precision tubing and High-Flo needle set should only be used to stop flow immediately in the case of an emergency. Use of the slide clamp may cause damage to the tubing and can affect the intended flow rate.
- Carefully inspect the Freedom60 Infusion Pump before use. Discontinue use of a pump that has been damaged, exposed to severe impact, or which fails to operate properly. Inspect the device for signs of wear (including cracks, fractures, corrosion, discoloration, pitting, or cracked seals) and discontinue use if device is unacceptably deteriorated.
- Avoid placing High-Flo needle set over a mole, tattoo, scar, muscle, or hardened or bruised areas, where proper needle set insertion could be difficult.
- Only perform infusions while stationary or walking. Infusions with movement other than walking may result in faster, slower, or more variable flow rates than specified. Testing has been performed to simulate walking and its effect on flow rates; the Freedom60 Infusion System has not been tested with any other physical activity.
- To obtain maximum accuracy of the pump, position the height of the pump within $\pm 7,6$ cm (3") of the infusion site, whether infusing in a stationary position or in motion. If the pump is higher than the infusion site, flow rate may increase. If the pump is positioned lower than the infusion site, flow rate may decrease.
- Ambient temperature may impact the flow rate. For most accurate results, it is recommended to infuse at a temperature of 20-25°C. Environmental pressure and humidity are not expected to significantly impact the flow rate.
- The shelf life of the device (prior to first use) has been verified as two (2) years after the manufacturing date when stored at 15-30°C at a humidity of 45-55%. Storage outside of these conditions, or use of the device after two (2) years of aging, may result in reduced performance or loss of device function. Refer to the device label for the manufacturing date of the device.

- The Freedom60 has been verified to withstand transport conditions of temperatures (-29-60°C) $\pm 2^{\circ}\text{C}$ and humidity of (15-85%) $\pm 5\%$ without compromising its functional performance. Transporting the device outside of these conditions may damage the device or result in unknown flow rate performance.

Freedom60® Diagram

Freedom60 Prefilled Syringe Adapter

Only for Hizentra® 50 ml prefilled syringes

Compatible Syringes

- **BD® Plastipak™ Luer-Lok® 50 ml**
EU Item Number: 300865 | US Item Number: 309653
- **Medline® 60 ml Luer Lock Syringe**
Item Number: SYR160010
- **Hizentra® 50 ml Single-Use Prefilled Syringe**
Item Number: 44206-0455-25

*Becton Dickinson & Co.

INSTRUCTIONS FOR USE

Before You Begin:

1. Gather your supplies and sanitize by cleaning your infusion work surface with antiseptic wipes or disinfecting solution. Wash your hands thoroughly and, if required, put on disposable gloves. Then lay out your supplies.
2. Make sure that the infusion pump's ON/OFF switch is in the OFF position and that the black tab within the syringe shield is at the end of its track. If the black tab is not touching the end of its track, fully turn the large wind knob clockwise.
3. Turn the pump ON and watch that the black tab moves smoothly along the full length of its track. Once complete, be sure to fully turn the large wind knob clockwise so the black tab is at the end of its track and you are ready for infusion.
4. Inspect the device to ensure all components are intact and free from cracks, damage, corrosion, or discoloration. Examine the inside of the syringe shield and ensure it is free of debris or contamination.
5. Verify that you are using the correct Precision tubing and High-Flo needle set prescribed by your healthcare provider. Inspect tubing and needle set packaging. Check that the bag is intact (no rips, tears, or holes) and fully sealed along all edges of the bag. If the packaging is damaged, do not use. Replace and contact your healthcare provider.
6. Ensure the medicinal product is at room temperature (20-25°C).

Before subcutaneous self-administration, patients and/or caregivers must be trained by a qualified healthcare provider.


Contact your healthcare provider if your pump is damaged, if the Precision tubing and/or High-Flo needle set packaging is damaged, or if you see any debris or contamination inside the pump.


If Using a Hizentra 50 ml Prefilled Syringe:

Do not use BD 50 ml or Medline 60 ml syringes with the Freedom60 Prefilled Syringe Adapter.


Installation Instructions for the Freedom60 Prefilled Syringe Adapter


Make sure adapter is in the correct orientation, with opening at the bottom.

Line the adapter up with the collar of the Freedom infusion pump.


Lightly squeeze the sides of the adapter and begin to push it into the collar.


Continue pushing the adapter into the collar until it clicks into place.

Removal Instructions for the Freedom60 Prefilled Syringe Adapter

Pinch the two retaining tabs inward by applying pressure to them.


While pinching tabs, pull the adapter out of the collar, opposite the direction it was inserted.

Step-by-Step Instructions for Subcutaneous Infusion

Part I – Prep System

1. Prepare Precision tubing & High-Flo needle set


Open the sterile package by tearing along the tear notch on the top of the package. Inspect needle set and tubing for damage. Check for bent needle set or damage to tubing. If damaged, replace and contact your healthcare provider.

2. Prepare syringe(s)


Refer to the medicinal product prescribing information or ask your healthcare provider for detailed filling instructions. Refer to **page 5** for a list of compatible syringes

- **Vial to syringe transfer [A]:** fill required dose into compatible syringe(s)
- **Prefilled syringe to syringe transfer:** fill required dose into compatible syringe(s) [B]
- **Hizentra® 50 ml prefilled syringe:** fits directly into pump with Freedom60 Prefilled Syringe Adapter, no transfers required. Refer to Freedom60 Prefilled Syringe Adapter Installation Instructions on **page 6**.

3. Attach tubing & needle set

Remove sterile caps from ends of the Precision tubing set and High-Flo needle set using care not to contaminate (touch) the ends. Connect the female end of the needle set to the male end of the tubing set [A]. Remove the cap from the female end of the tubing and attach the female end with luer disc to the syringe [B].

4. Prime (fill) Precision tubing

Always follow your healthcare provider's protocol.

NOTES:



- It is recommended to insert the needle set dry, meaning no medicinal product is present on the needle tip(s), to minimize infusion site irritation.
- To best see the medicinal product in the tubing during priming, it is recommended to prime the tubing against a dark, solid-colored surface in a well-lit area.

There are two options for priming: by hand or by pump. See step by step instructions below.

Priming by hand:


Lay the tubing and needles on work surface.

- Push the syringe plunger and follow the medicinal product as it flows through the tube.
- When the fluid reaches the desired point, release pressure from the plunger to stop the flow.
- Focus on a single needle and try to stop the flow when the fluid approaches the needle. Be careful not to prime to the needle tip.


Priming by pump:

Make sure the pump is in the OFF position and the black tab inside the clear syringe shield is at the end of its track. If the black tab is not at the end of its track, turn the large wind knob clockwise. *You should not need to use significant force to load or remove the syringe. If you are having trouble, stop the pump and make sure the black tab is at the end of its track before removing the syringe.*

- If using a Hizentra 50 ml prefilled syringe, be sure the Freedom60 Prefilled Syringe Adapter is installed prior to loading syringe into pump as directed below (see **page 6** for instructions). If not using a Hizentra 50ml prefilled syringe please be sure the Freedom60 Prefilled Syringe Adapter is removed before moving forward.
- With syringe gradations facing up, load the assembled syringe into the pump and ensure that the luer disc is fully seated in the pump nose.
- Lay Precision tubing and High-Flo needles on work surface. Turn the pump ON to prime (fill) the tubing and watch the tubing as the medicinal product approaches the needle. When the medicinal product reaches the desired point, turn the ON/OFF switch to the OFF position and immediately turn the large wind knob clockwise to release pressure on the plunger.


Part II – Prep Sites

5. Prepare sites

- Select, clean and let site(s) dry before inserting needle set. Carefully remove the shield from the needle tip, with care not to touch the needle.
- Always refer to the medicinal product prescribing information and recommendations from your healthcare provider for infusion site location(s). Common areas for subcutaneous infusion include the abdomen, thighs, side of the upper hips and back of the arms.


*Avoid placing needle over a mole, tattoo, scar, muscle, hardened or bruised areas, where proper needle insertion could be difficult.


6. Insert needles

- Pinch the skin and insert each needle into the subcutaneous tissue at a 90° angle. Secure the needle (following the instructions in Step 7) before inserting a needle in another site.


7. Secure needle set

- Peel the printed side from the dressing to expose adhesive.
- Secure the needle by placing the adhesive dressing in the center of the needle butterfly. Smooth it outward over skin.
- Repeat for each infusion site.

8. Check for blood return

- Check for blood return if instructed by your healthcare provider by gently pulling back on the syringe plunger **[A]**. Watch to make sure no red/pink appears in tubing near your sites.
- If blood return exists and if instructed by your healthcare provider, either clamp the flow to the needle site(s) **[B]** or remove all needles, attach a new needle set, and start again from Step 3.

Part III – Infuse

When using the Hizentra 50 ml prefilled syringe, be sure the Freedom60 Prefilled Syringe Adapter is installed prior to inserting the syringe into the pump.

When NOT using a Hizentra 50 ml prefilled syringe, be sure to remove the Freedom60 Prefilled Syringe Adapter. **See page 6 for instructions.**


9. Insert syringe & turn ON

- Ensure the black tab is at the back of the pump. With syringe gradations facing up, insert the syringe into the pump.
- Ensure the luer disc is fully seated in the pump's nose.
- Turn pump ON

Note: To reduce flow rate variability, try to keep the pump level with your infusion sites.


Caution: if luer disc is not fully seated in pump nose it may eject from the pump.

10. Infuse

- Periodically check that the pump is working properly by seeing that the syringe plunger is moving.

If using multiple syringes:

- Once the first syringe is empty, turn the pump OFF and wind the black tab to the end of its track. Remove the syringe from the pump and disconnect from tubing. Using care not to contaminate the ends, connect the additional syringe to the luer disc end of the Precision tubing.
- When using the Hizentra 50 ml prefilled syringe, be sure the Freedom60 Prefilled Syringe Adapter is installed prior to inserting the syringe into the pump (see **page 6** for instructions).
- When NOT using a Hizentra 50 ml prefilled syringe, be sure to remove the Freedom60 Prefilled Syringe Adapter (see **page 6** for instructions).
- Load the prepared syringe into the pump. Turn the pump ON to continue infusion. Repeat until total dosage is complete.

If you are experiencing severe pain or discomfort, stop use immediately and contact your healthcare provider.

To pause infusion immediately, turn the pump OFF and wind the large knob until the black tab is at the end of its track.


If the device malfunctions or is not performing as expected, refer to the Troubleshooting instructions on **page 12**. If problems persist, discontinue use and consult your healthcare provider.

Any serious incident should be reported to your healthcare provider and KORU Medical Systems, Inc. at +1-800-624-9600. For patients and/or users in the European Union (EU): Any serious incident should be reported to the competent authority in your region.

Part IV – Complete Infusion

11. Turn off & wind back

- When the syringe is completely empty and total dose has been administered, turn the pump OFF. Wind the large knob until the black tab is at the end of its track.


12. Remove syringe

- Remove empty syringe and its tubing.

13. Remove High-Flo Needle & cleanse sites

- Holding the needle in place, peel back the surrounding adhesive dressing and remove.
- Pull needle straight out to remove. To use safety feature, close wings over the needle and snap shut.
- Cleanse each site and cover with a bandage if needed.


14. Discard disposables

Discard your supplies in a secure sharps container according to your healthcare provider's instructions.

15. Clean and disinfect device

- Clean and disinfect the pump as soon as possible after use and avoid delays between cleaning/ disinfection steps. See **page 14** for full cleaning and disinfection instructions.

TROUBLESHOOTING

If the suggestions in this section do not solve your problem, or if problems persist, discontinue use and consult your healthcare provider. Any serious incident should be reported to your healthcare provider and KORU Medical Systems, Inc. at **+1-800-624-9600**. For patients and/or users in the EU: Any serious incident should be reported to the competent authority in your region.

If you are experiencing pain or extreme discomfort, stop the infusion immediately and call your healthcare provider.

Syringe will not load or remove from Infusion Pump:

NOTE: You should not need to use significant force to load or remove a syringe.

- Make sure the pump is in the OFF position and that the black tab is at the end of its track. If the black tab is not at the end of its track, fully wind the large knob clockwise and try removing syringe again.
- Make sure you are using a syringe listed on **page 2**. These are the only syringes that have been tested for use with the Freedom60.
- If using a Hizentra® 50 ml prefilled syringe, ensure the Freedom60 Prefilled Syringe Adapter is installed prior to loading the syringe (see **page 6** for instructions)
- When NOT using a Hizentra 50 ml prefilled syringe, be sure to remove the Freedom60 Prefilled Syringe Adapter (see **page 6** for instructions)

Syringe will not stay inside in the Infusion Pump:

- Make sure you are using a Precision Flow Rate Tubing set and that the luer disc end of the tubing has been connected to a compatible syringe (see **page 2** for details).
- Make sure you are using a syringe listed on **page 2**. These are the only syringes that have been tested for use with the Freedom60.
- Make sure the flange shape of the syringe is seated securely into the pump's nose.
- Make sure the luer disc on the tubing set is seated properly in the nose of the pump.
- Make sure you have attached the syringe directly to the tubing and NOT directly to the subcutaneous needle set.

Medicinal product (5 ml or less) is left in the syringe:

- Make sure you are using a syringe listed on **page 4**. These are the only syringes that have been tested for use with the Freedom60. Contact your healthcare provider if the syringe does not completely empty.

Medicinal product is not flowing through the tubing:

- Make sure that the pump is in the ON position.
- Make sure all the tubing is not clamped.
- Using care not to touch or contaminate the ends, disconnect the tubing set from the needle set, and check for medicinal product drip.
- If the medicinal product does not drip, replace the tubing as it may be damaged.

The flow is slow (infusion is taking longer than expected):

- Flow rate can vary for each patient and each infusion. Administration may be slow based on how well the medicinal product is absorbed through the tissue. Some infusions may be faster than others. The first infusions may take longer than expected because the body may need to adapt.
- Verify that the tubing is not kinked, damaged, or punctured. If the slide clamp has been used, the tubing may be damaged and should be replaced.
- Make sure you are using a syringe listed on **page 2**. These are the only syringes that have been tested for use with the Freedom60.
- Ensure the pump is level with the infusion sites. If the pump is lower than the infusion sites the flow rate may be slower.

- Only perform infusions while stationary or walking. Infusions with movement other than walking may result in faster, slower, or more variable flow rates than specified.
- Ensure that your medicinal product is at room temperature and that the ambient temperature in your infusion space is between 20-25°C. Infusing at a lower temperature will result in a slower flow rate.
- Avoid placing needle set over a mole, tattoo, scar tissue, muscle, or hardened or bruised areas. The tissue may not absorb the fluid as intended.
- If you prefer a faster flow, you may need more sites, a different needle set and/or a faster flow rate tubing set. Contact your healthcare provider.
- Your pump may have reached the end of its use life. The spring mechanism in the pump may be worn out, resulting in a slower flow than expected. The test procedure on **page 15** can be performed to verify if the device is performing according to its specifications. Contact your healthcare provider for more information on replacing your infusion pump.

The flow is faster than expected (infusion is taking less time than expected):

- Ensure the pump is level with the infusion sites. If the pump is above than the infusion sites the flow rate may be faster.
- Make sure you are using a syringe listed on **page 2**. These are the only syringes that have been tested for use with the Freedom60.
- Flow rate can vary for each patient and each infusion. Administration may be faster based on how well the fluid is absorbed through the tissue. Some infusions may be faster than others. Subsequent infusions after initial dosing may flow faster as your body adapts to the infusions.
- Ensure you have connected your flow rate tubing. If the syringe is connected directly to the needle set, there will be no control over the flow rate, and the syringe may eject from the pump.
- Only perform infusions while stationary or walking. Infusions with movement other than walking may result in faster, slower, or more variable flow rates than specified.
- Ensure that your medicinal product is at room temperature and that the ambient temperature in your infusion space is between 20-25°C. Infusing at a higher temperature may result in a higher flow rate.
- If you prefer a slower flow, you may need fewer sites, a different needle set and/or a slower flow rate tubing set. Contact your healthcare provider.

Flow is continuing after infusion has stopped:

- This is a normal function of the pump as it is designed to maintain pressure during and after the infusion to prevent medicinal product backflow.
- To release pressure from the syringe plunger and to stop the flow, turn the large wind knob clockwise so that the black tab is at the end of its track.
- You can also use the slide clamp to stop the flow immediately. This should only be done in the case of an emergency as use of the slide clamp can damage the tubing. After using the slide clamp, replace the tubing set and/or needle set before continuing the infusion.

There are local site reactions (e.g., swelling, redness or soreness):

- Allowing medicinal product to contact the needle when priming may cause irritation. It is recommended to insert subcutaneous needle set dry as the medicinal product may irritate the skin.
- Make sure that the needle set(s) are long enough to reach the subcutaneous layer. If the selected needle is too short, leaking at the site may occur. Contact your healthcare provider about using a longer needle.
- Make sure that the needle set(s) are not too long, as they may hit muscle. Contact your healthcare provider about using a shorter needle.
- Ask your healthcare provider about trying a slower flow rate tubing set as the rate may be too fast.
- Rotate infusion sites if recommended by your healthcare provider. Periodically returning to sites that worked well in the past may provide best results.

CARE, MAINTENANCE AND DISPOSAL:

The Freedom60® Infusion Pump does not require any calibration to be performed by the user.

Transport:

The Freedom60 has been verified to withstand transport conditions of temperatures (-29-60°C) $\pm 2^{\circ}\text{C}$ and humidity of (15-85%) $\pm 5\%$ without compromising its functional performance. Transporting the device outside of these conditions may damage the device or result in unknown or inaccurate flow rate performance.

Storage:

The Freedom60 Infusion Pump has been verified to maintain its performance when stored for up to two (2) years at a temperature of 20-25°C and humidity of 45-55%. Storing the device at a higher or lower temperature or in a more humid environment may damage the pump and/or reduce the expected lifetime.

Expected Lifetime and Shelf Life:

The Prefilled Syringe Adapter has been verified for a minimum use life of 360 cycles (15 uses per month for two years). Inspect the device prior to use. Do not use if the Adapter is cracked, deformed, worn down, or otherwise damaged. Discard and contact your healthcare provider for a replacement if the Adapter is not usable.

The Freedom60 Infusion Pump has an expected shelf life of two (2) years prior to first use, and has been tested to meet its stated performance specifications for 360 uses (representative of two (2) years of use life when used 15 times per month) when stored and cleaned according to the instructions below.

Inspect the device prior to use and after each use, and monitor its performance during infusion. Contact your healthcare provider for a replacement if the device does not perform as intended or if the device shows signs of deterioration such as corrosion, discoloration, pitting, and cracks.

Disposal Instructions for Patients:

Contact your healthcare provider for specific disposal instructions for the infusion pump for your region.

Discard used infusion accessories (needle sets, tubing sets, and other infusion supplies) in a sharps container according to your healthcare provider's instructions. Ensure the needle set butterfly closure is securely fastened over the needle before disposal.

Cleaning and Disinfection Instructions:

Between uses, the infusion pump must be cleaned and then disinfected.

After cleaning and disinfection, inspect the device for unacceptable deterioration such as corrosion, discoloration, pitting, and cracks. Monitor your device for unexpected or deteriorating performance (e.g., slower than expected flow rate or stoppage of flow). If your device is damaged or malfunctioning, contact your healthcare provider for instructions on replacement and disposal.

Caution: Clean only those areas that are exposed and external. No attempt should be made to clean any part of the infusion pump that is not easily accessible.

Cleaning Procedure:


1. Clean the pump with a soft cloth dampened with a weak mixture of mild detergent and warm water (minimum ratio of one (1) part detergent to 50 parts water by volume).
2. Using the prepared detergent solution and a clean non-linting wipe or soft cloth, wipe all the external surfaces of the pump, including the pump nose and syringe tray up to the syringe shield for at least one (1) minute. Ensure the device is fully open to allow all surfaces to be cleaned. During the one (1) minute wipe, pay special attention to the ridges, crevices, raised lettering during wiping. Replace soiled cloths or wipes as needed, changing wipes when necessary to ensure that all surfaces are cleaned.

3. Using a clean non-linting wipe or soft cloth wetted with room temperature tap water (wet but not dripping), wipe all the external surfaces of the pump, including the pump nose and syringe tray up to the syringe shield. Pay special attention to the ridges, crevices, raised lettering during wiping. Continue wiping until all residue is removed to ensure the pump is thoroughly clean. Replace or re-wet cloth or wipes as needed, changing wipes when necessary to ensure that all surfaces are rinsed.
4. Dry the pump using a clean non-linting wipe or soft cloth.
5. Inspect the pump for any visible residue and debris to ensure that the device is thoroughly cleaned. If the device has remaining visible residue and debris following cleaning, repeat the cleaning steps (1 through 4). Once all visible residue and debris is removed, disinfect pump per the disinfection procedure.

Disinfection Procedure:

1. Use pre-saturated IPA wipes, or a non-linting cloth or wipe saturated with 70% Isopropyl Alcohol (IPA) (wetted but not dripping) to thoroughly wipe all exterior surfaces of the Freedom60 Infusion Pump.
 - Ensure all external surfaces of the pump and adapter, including the pump nose and syringe tray up to the syringe shield, are wiped. Pay special attention to the ridges, crevices, raised lettering during wiping. Ensure the device is fully open to allow all surfaces to be cleaned.
-  **Caution:** Clean only those areas that are exposed and external. No attempt should be made to clean any part of the pump that is not easily accessible.
- Allow all surfaces to remain visibly wet for a minimum of five (5) minutes. During the five (5) minute contact time, use additional wipes to ensure all contacted surfaces remain wet for the full contact duration time.
2. Thoroughly dry the pump using non-linting wipe(s) or allow to air dry.
 3. Visually inspect the pump for signs of damage or wear (such as corrosion, discoloration, pitting, or cracks.)

DEFINITION OF SYMBOLS

	Catalog Number		Consult IFU before use. IFU available electronically at the indicated location.
	Manufacturer		For Prescription Use Only
	Authorized Representative in the European Community		Device is not safe for use in a magnetic radiation environment
	Swiss Representative		Do Not Reuse
	UK Responsible Person		Do Not Re-sterilize
	Lot Number		Item is a medical device
	Serial Number		Single sterile barrier packaging system
	Unique Device Identification (GTIN)		Do not use if package is damaged
	Expiration Date		CE Mark Notified Body Number (BSI NL: 2797)
	Manufacturing Date		Not manufactured with natural rubber/latex
	Temperature Range to which device can be exposed		Importer
	Humidity Range to which device can be exposed		Caution required when using device. Refer to IFU
	Sterilized with Gamma Radiation		

For warranty information, please visit www.korumedical.com/warranty and contact your healthcare provider.

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This page was intentionally left blank

This page was intentionally left blank

KORU Medical Systems, Inc.
100 Corporate Drive, Mahwah,
NJ 07430 USA.
+1 (845) 469 2042
800-624-9600
korumedical.com

ICON (LR) Limited
South County Business Park,
Leopardstown, Dublin 18,
D18 X5R3, Ireland

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

MedEnvoy
Prinses Margrietplantsoen 33
Suite 123
2595 AM The Hague
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

1-800-624-9600 | korumedical.com | @korumedical

The Freedom60® Syringe Infusion System, Precision™ Flow Rate Tubing, Precision™ Flow Rate Controller, and High-Flo™ SubQ Needle Sets are trademarks of KORU Medical Systems, Inc. All other trademarks are the property of their respective owners.
©2025 KORU Medical Systems, Inc.; All Rights Reserved.

348442-B Rev A - 08/2025